

CYNOSURE®

Laser Lipolysis Redefined

Smartlipo[™] Triplex



Only Smartlipo TriPlex offers 3 wavelengths

Including 1440 nm for faster disruption of adipose tissue – delivered by multiplex sequencing and intelligent delivery systems. It's laser lipolysis redefined, only from Cynosure. Why settle for anything less?

Greater absorption.

1440 nm has 40 times greater absorption in fatty tissue than diode lasers.





More choice.

Multiplexed combination wavelength therapy for optimal fat disruption, high-definition sculpting capability and tissue tightening through coagulation.

Better Results.

SmartSense, ThermaGuide and ThermaView intelligent delivery systems help you treat with the maximum dose, safely.





Dramatic
Results in
Less Time.

No wonder
Cynosure's Smartlipo
multiplexing
technology is the
leading choice for
laser lipolysis around
the world.

Smartlipo TriPlex™

New Improved Design!

Longer Lasting
ThermaGuide Design

Easier to Use
Software Design

Larger Display

Choose the right energy level and wavelength without overtreating or under treating.

For the most accurate and efficient laser lipolysis therapy available today choose Smartlipo TriPlex.

Patient's growing desire for both a youthful appearance and minimally invasive procedures will drive the laser assisted lipolysis market through the next 5 years¹.

Be ready to meet your patient's ever-increasing demands with the expanded clinical capabilities of Cynosure's Smartlipo TriPlex – the first multiplex triple-wavelength laser lipolysis workstation that's designed to treat all your patient's body contouring concerns:

- More effective laser lipolysis
- Proven tissue coagulation resulting in tightening
- Reduced treatment time
- Less Downtime
- High-definition sculpting capabilities



Be Smart. Be Sure.
Pick the one that's just right for you.

With more wavelengths in a single workstation, you can treat a wider variety of patients more effectively than ever before. The Smartlipo TriPlex all-in-one laser system allows you to:

Personalize your treatment for:

- **More efficient LaserBodySculpting™ for faster treatments**
- **High definition treatments**
- **Targeted heating for faster recovery**
- **Areas of skin laxity**
- **Broad tissue heating for proven tissue tightening results through coagulation**
- **Smaller treatment areas, like the face and neck**
- **Treat axillary sweat glands associated with hyperhidrosis**
- **Easily address more fibrous tissues associated with lipomas and gyneconastia**



“The biggest advantage to Smartlipo TriPlex is speed. If I can do more cases or complete them faster, I can reduce my overhead, treat more areas at one time and be more profitable.”

-John Decorato, MD
Staten Island, NY

Proven to Drive Patients to Your Practice.



Patient demand for Smartlipo is unprecedented and growing more rapidly than ever. Smartlipo's superior satisfaction rate – 94% of patients treated would recommend Smartlipo² - has made Smartlipo the most searched for cosmetic procedure on the internet.

- Smartlipo Patient Website brings patients to you. www.smartlipocynosure.com
- Exclusive advertiser on www.liposuction.com
- Smartlipo was voted best overall procedure by patients in 2009 and 2010³.
- 82% of all laser lipolysis internet research leads patients to Smartlipo³

Realize an immediate return on your Smartlipo TriPlex investment by capitalizing on all the marketing opportunities Cynosure has to offer.



“ My experience with Smartlipo has been very positive. Without a doubt, I have seen less bruising, less swelling, and faster recovery than liposuction alone. Smartlipo also has a strong market value. 125 patients in the last 15 months have been to my practice specifically seeking Smartlipo.”


- Joseph Hunstad, MD
Huntersville, N.C.



Better Science. Better Results.

Smartlipo TriPlex claims are supported by more than 9 years of clinical studies and science. This substantial body of evidence provides you with a scientific basis to merge Smartlipo TriPlex with your techniques for the most ideal treatment outcomes.

To get your free copy of a compilation of Smartlipo clinical evidence, please visit www.cynosure.com/proof



Intelligent delivery systems help you deliver the maximum dose, safely.

SmartSense. ThermoGuide. ThermoView. Cynosure's unique, intelligent delivery systems are so smart, they can help you prevent over-treatment, determine your endpoint and minimize complications with integrated temperature and motion sensing.

You'll find all three systems are available for Smartlipo products, the laser lipolysis system that started it all.

| | | | |
|-------------------------|---|-------------------|--------|
| System Type | Nd:YAG | | |
| Wavelength | 1064nm | 1320nm | 1440nm |
| Pulse Width (Max) | 500µs | 1000µs | 1000µs |
| Power Output Options | 40W | 24W | 15W |
| Max Rep Rate | 40Hz | 40Hz | 25Hz |
| | MultiPlexing Options | | |
| | 1064nm and 1320nm | 1064nm and 1440nm | |
| | 500µs and 1000µs | 500µs and 1000µs | |
| | 40Hz | 40 Hz and 25Hz | |
| | | | |
| Aming Beam | Red | | |
| Electrical Requirements | 230 VAC/16A | | |
| Size (HxWxD) | 37.4" x 19.7" x 25.6" (95cm x 50cm x 65cm) | | |
| Weight | 220 lbs (100kg) | | |

Please contact your Smartlipo TriPlex Surgical Laser Consultant at 800.886.2966 or visit www.cynosure.com.

Sources:

1. Millenium Research Group Report "North America Markets for Body Contouring Devices 2009"
2. www.americanhealthandbeauty.com
3. According to www.americanhealthandbeauty.com

Before and After Photos:

1. Courtesy of C. Petti, MD
2. Courtesy of J. Millard, MD
3. Courtesy of T. Simon, MD
4. Courtesy of B. DiBernardo, MD



www.cynosure.com

CORPORATE HEADQUARTERS

United States
 Cynosure, Inc.
 5 Carlisle Road
 Westford, MA 01886 USA
 Tel: 978.256.4200
 Toll-Free: 800.886.2966

